


2x2 Diamond 38.5" x 38.5"
(Template at full size)

SKIN A


Indicates 1.5" bleed
Indicates edge of skin
Indicates 3" focus area




1x1 Flat Front 28" x 28"
(Template at full size)

SKIN B

Indicates 1.5" bleed
Indicates edge of skin
Indicates 3" focus area




Will be applied
to frame like this




1x1 Flat Front 28" x 28"
(Template at full size)

SKIN C

Indicates 1.5" bleed
Indicates edge of skin
Indicates 3" focus area




Will be applied
to frame like this




1x1 Single Twist 28.5" x 28.5"
(Template at full size)

SKIN D

Indicates 1.5" bleed
Indicates edge of skin
Indicates 3" focus area




Will be applied
to frame like this




1x1 Flat Front 28" x 28"
(Template at full size)

SKIN E

Indicates 1.5" bleed
Indicates edge of skin
Indicates 3" focus area




Will be applied
to frame like this




1x1 Flat Front 28" x 28"
(Template at full size)

SKIN F

Indicates 1.5" bleed
Indicates edge of skin
Indicates 3" focus area




Will be applied
to frame like this




1x1 Flat Front 28" x 28"
(Template at full size)

SKIN G

Indicates 1.5" bleed
Indicates edge of skin
Indicates 3" focus area




Will be applied
to frame like this



1x1 Flat Front 28" x 28"
(Template at full size)

SKIN H

Indicates 1.5" bleed
Indicates edge of skin
Indicates 3" focus area



Will be applied
to frame like this